



VBM30-04-0041

ADULT QUICKTRACH

- Urgent CONIOTOMY sterile set QUICKTRACH type
- Complete with:
 - o 1 plastic canula with intubation internal-blade stylet with diameter of 4 mm
 - o 1 ISO joint for the connection to resuscitation bag
 - 1 sharp conical-tip intubation stylet that avoids windpipe incision using classic scalpel and a "lock" that avoids penetrating too much in depth
 - o 2 lateral wings allowing airways safe fixing

DIR/QM - Rev. 02_11/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.