



**MA2023001**

## **INTUBATION STYLET 2.0 MM**

- Tracheal intubation stylet made of flexible, plastic-coated aluminium, single-use and sterile.
- Rounded distal tip to reduce damage to the tracheal mucosa.
- Non-conductive, latex- and phthalate-free.
- Dimensions: D.E. 2.0 mm. Length 272 mm. D.I. 2.5 ÷ 4.0 mm.
- EC medical device.

*DIR/QM – Rev. 02\_06/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.