



**MA3034001**

### **INTUBATION STYLET 3.3 MM**

- Tracheal intubation stylet made of flexible, plastic-coated aluminium, single-use and sterile.
- Rounded distal tip to reduce damage to the tracheal mucosa.
- Non-conductive, latex- and phthalate-free.
- Dimensions: D.E. 3.3 mm. Length 346 mm. D.I. 4.5 ÷ 5.5 mm.
- EC medical device.

*DIR/QM – Rev. 02\_06/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.