



MA5040001

INTUBATION STYLET 4.7 MM

- Tracheal intubation stylet made of flexible, plastic-coated aluminium, single-use and sterile.
- Rounded distal tip to reduce damage to the tracheal mucosa.
- Non-conductive, latex- and phthalate-free.
- Dimensions: D.E. 4.7 mm. Length 405 mm. D.I. 6,0 ÷ 10,0 mm.
- EC medical device.

DIR/QM – Rev. 02_06/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.