



ETAN30C02

3.0MM CUFFED ORO-NASAL TRACHEAL TUBE

- Disposable, cuffed, oro-nasal tracheal tube
- Made of thermosensitive polyvinyl chloride (PVC)
- Intended for the following uses:
- positive pressure ventilation;
- induction of anaesthesia if it is not possible to use a face mask;
- removal/protection of secretions from the trachea and bronchi during suction.
- 15M connector
- Graduated with depth marking to check positioning in relation to the vocal cords
- Radiopaque line for clear identification under X-rays
- Equipped with:
- Low pressure, high volume, atraumatic cuff
- Blue indicator balloon with Luer fitting and one-way valve
- Atraumatic, rounded edge Murphy eye
- Atraumatic, pre-curved Magill tip in compliance with applicable ISO/EN standards.

DIR/QM - Rev. 01_31/03/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.