



ETAN35C02

3.5MM CUFFED ORO-NASAL TRACHEAL TUBE

- Disposable, cuffed, oro-nasal tracheal tube
- Made of thermosensitive polyvinyl chloride (PVC)
- Intended for the following uses:
 - positive pressure ventilation;
 - o induction of anaesthesia if it is not possible to use a face mask;
 - o removal/protection of secretions from the trachea and bronchi during suction.
- 15M connector
- Graduated with depth marking to check positioning in relation to the vocal cords
- Radiopaque line for clear identification under X-rays
- Equipped with:
 - o Low pressure, high volume, atraumatic cuff
 - o Blue indicator balloon with Luer fitting and one-way valve
 - o Atraumatic, rounded edge Murphy eye
 - o Atraumatic, pre-curved Magill tip in compliance with applicable ISO/EN standards.

DIR/QM - Rev. 01_31/03/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.