



**7 mm**

**NC7015001**

### **NASOPHARYNGEAL CANNULA I.D. 7 MM**

- Device for maintaining patency of the nasopharyngeal cavity
- For nasopharyngeal intubation during anesthesia, intensive care, and emergency situations.
- Made of PVC, soft, flexible, resistant to collapse
- Rounded edges, smooth atraumatic surface
- 15M connector for connection to ventilation systems
- Supplied individually in sterile packaging
- Latex-free and phthalate-free (DEHP/DOP free)
- Single use
- Internal diameter: 7.0 mm
- Length: 150 mm

*DIR/QM – Rev. 02\_12/03/2026*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.