



#### FT100B

#### INFRARED DIGITAL THERMOMETER\*

- LCD display with 3½ figures
- Temperature unit of measurement: Celsius degrees (°C) and Fahrenheit (°F)
- Measurement range: 32.0 – 42.9 °C (89.6 – 109.2 °F)
- Saves last 10 performed measurements
- Powered by not rechargeable battery, PC 2032 type (1pc)
- The batteries last about 2500 uses
- Size: 105x35x35mm
- Net weight: 36gr
- Acoustic alarm: beep sound for ON-OFF
- Lifetime: 3 years
- Complete with protective plastic case

DIR/QM – Rev. 02\_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.