



**OXI0587**

## **FLOWMETER 0-15 LPM AFNOR PLUG MD**

- Flowmeter for oxygen supply
- Output flow between 0-15 lpm (liters per minute)
- Real-time indication of oxygen supply through graduated scale placed on the burette
- Nickel-plated brass body
- AFNOR plug for direct connection to the terminal unit
- 9/16" connection to the humidifier or to the hose connector (optional)
- Filter positioned at the inlet of the burette to protect the patient and the flowmeter itself
- Gas flow adjustment is possible through the knob on the front of the device
- Compatible with disposable or reusable humidifiers
- Made in Italy
- Class IIa medical device in accordance with Directive 93/42/EEC

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.