



**OS/50P**

## **PAEDIATRIC OXYGEN THERAPY MASK WITH RESERVOIR**

- Oxygen therapy mask with reservoir and anti-backflow valve
- High concentration (60% Oxygen)
- Manufactured in transparent material, atoxic and soft, complete with adjustable elastic
- Metal nose clip allowing a safe and comfortable appliance
- Supplied with 210 cm crushproof tube and universal connector
- Size: pediatric

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.