



OXIP7010P

OXYGEN CYLINDER EXCHANGER WITH PLATE

- High performance mechanical exchanger with lateral or rear inlets/outlets (prepared by the manufacturer)
- Accommodates the pressure gauge inside the exchange knob for pressures up to 15 bar (1500 kPa)
- Thanks to the special construction, the pressure gauge is protected against accidental mechanical shocks and impacts
- The pressure gauge fitted inside the exchanger allows the size of the oxygen bar to be reduced
- 3/8" inlet and outlet connections perfectly compatible with oxygen hoses manufactured by Boscarol
- Allows a pressure transducer to be fitted on the outlet for measuring the nominal pressure of the kit
- The exchanger includes amounting plate for flush mounting in the ambulance.
- The plate is made of stainless steel and can be painted according to customer requirements
- The Boscarol cylinder exchanger is part of the oxygen kit and should be considered as a spare part

DIR/QM – Rev. 01_02/01/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.