



## BSU402

### OB500 MOD FA MSU

- OB500 medical suction unit equipped with autoclavable collection jar OB-J FA 1000 ml with overflow valve and protective filter inserted directly in the lid, autoclavable max 30 cycles
- The device is designed for medical vehicles and is stationary, not portable
- Nominal suction power: 800 mbar (80 kPa)  $\pm 10\%$
- Nominal suction flow: 30 LPM at free air speed  $\pm 10\%$
- Maximum operating time (free-cycle): 60 minutes  $\pm 10\%$
- Degree of protection against infiltration of liquids and solids: IP32d
- Functioning: 12÷15 Vdc
- Maintenance-free piston pump
- Vacuum control and regulation unit with ABS module
- Dimensions of suction unit: 175x175x100 mm
- Dimensions of ABS module adjustment unit: 90x130x85 mm
- Weight: max. 2 kg including mounting parts
- Class IIa medical device in compliance with Regulation (EU) 2017/745 (MDR), all applicable regulations and main reference standards

DIR/QM – Rev. 02\_12/06/2023

User manual available on the [eIFU](#) webpage

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.