



**BSU506**

### **OB-J JAR 1000ML FOR DISPOSABLE BAGS**

- Jar OB-J for secretion collection made of polypropylene for medical use
- Consisting of a jar and an adapter that allows the jar to be used with Serres' LINER disposable bags (not included)
- Complete with "L" joint
- The canister is vapor autoclavable at a max temperature of 121°C and at a pressure of 2 bar (200kPa)
- The jar can also be disinfectable accordingly to local hospital procedures

*DIR/QM – Rev. 02\_12/06/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.