



## BSU810

### WALL BRACKET OB WB

- Wall bracket anchorage for OB2012, OB3000 and OB1000 medical suction units
- Integrated device charging system
- Anchorage system tested to 10 g in accordance with the current EN 1789 standard
- Extremely compact and functional holder not protruding from the wall of the rescue vehicle, with quick release button
- 12Vdc power supply for OB2012 and OB1000; 11÷30Vdc for OB3000
- Weight: about 0.8 kg
- Dimensions: about 240 x 255 x 5 mm

DIR/QM – Rev. 03\_22/03/2024

User manual available on the [eIFU](#) webpage

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.