



#### M03.1.003

### PROTECTION FILTER FOR SUCTION UNIT\*

- Disposable protection filter with PTFE filter holder and hermetic styrene container
- Maximum liquid inlet pressure (tested with water): 16 psi
- Filtering and antibacterial capacity: >99.9999%. Filter area: 24.9 cm<sup>2</sup>.
- Device packaged and sold individually.
- Can be used as an alternative to the BSU999 filter as also stated in the IFU of the device
- EC Medical device

DIR/QM - Rev. 02\_27/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.