



**SC125003**

## **SUCTION UNIT CATHETER CH12 WHITE**

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
  - Length: 50 cm.
  - Size: CH 12.
  - Outer diameter: 4.0 mm.
- Colour: white.
- EC Medical Device

*DIR/QM – Rev. 02\_27/07/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.