



SC145003

SUCTION UNIT CATHETER CH14 GREEN

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
 - Length: 50 cm.
 - Size: CH 14.
 - Outer diameter: 4.7 mm.
- Colour: green.
- EC Medical Device

DIR/QM – Rev. 02_27/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.