





## SC165003

## **SUCTION UNIT CATHETER CH16 ORANGE**

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
  - $\circ\,$  Length: 50 cm.
  - o Size: CH 16.
  - o Outer diameter: 5.3 mm.
- Colour: orange.EC Medical Device

DIR/QM - Rev. 02\_27/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.