





SC205003

SUCTION UNIT CATHETER CH20 YELLOW

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
 - $\circ\,$ Length: 50 cm.
 - o Size: CH 20.
 - o Outer diameter: 6.7 mm.
- Colour: yellow.
- EC Medical Device

DIR/QM - Rev. 02_27/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.