



**SC105003**

## **SUCTION UNIT CATHETER CH10 BLACK**

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
  - Length: 50 cm.
  - Size: CH 10.
  - Outer diameter: 3.3 mm.
- Colour: black.
- EC Medical Device

DIR/QM – Rev. 02\_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.