



FS-6540

HAND SUCTION UNIT

- 300 ml jar for secretion collection
- Portable, functional and safe suction unit, operated exclusively by hand pressure, allowing to achieve a depression up to 450 mbar
- Its "GUN" shape permits to use just one hand during emergency activities
- Dimensions: approx. 22x8x18cm
- Weight: approx. 275 g
- It is a CE medical device

DIR/QM – Rev. 01_20/01/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.