



VAL40902

RESUSCITATION CASE SILICONE BAG

- ABS resuscitation case with wall anchorage containing:
 - 1 autoclavable silicone resuscitator bag
 - 1 Reservoir and 1 oxygen tube
 - 3 silicone transparent ventilation masks
 - 5 Guedel airways
- Dimensions: about 39,5x13,5x26(h) cm
- Weight: about 2,3 Kg

DIR/QM – Rev. 02_30/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.