



VAL40910

REANIMATION SET WITH SILICONE BAG

- Reanimation bag made of waterproof, tear-proof fabric, with carrying strap, featuring:
 - transparent autoclavable silicone resuscitator bag with turning connection security valve
 - 2 transparent ventilation masks with silicone board
 - reservoir
 - connecting tube for oxygen
 - 5 Guedel airways
- Dimensions: about 28x15x16 cm
- Weight: about 1 kg

DIR/QM – Rev. 02_30/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.