



VAL40916

REANIMATION SET WITH PVC BAG

- Reanimation Bag made of waterproof, tear-proof fabric, with carrying strap featuring:
 - disposable resuscitator bag with turning connection security valve
 - 2 transparent ventilation masks
 - reservoir
 - connecting tube for oxygen
 - 5 Guedel airways
- Dimensions: about 28x15x15h cm
- Weight: about 1 kg

DIR/QM – Rev. 02_30/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.