



**VAL40905**

## **RESUSCITATION CASE PVC BAG**

- ABS resuscitation case with wall anchorage system containing:
  - 1 pvc disposable resuscitator bag
  - 1 reservoir and 1 oxygen tube
  - 3 disposable ventilation masks
  - 5 Guedel airways
- Dimensions: about 39,5x13,5x26(h) cm
- Weight: about 2,2 kg

*DIR/QM – Rev. 02\_30/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.