



DIA02009

ANEROID TABLE SPHYGMOMANOMETER*

- Circular wall and/or table aneroid sphygmomanometer (diameter 170mm)
- Scale from 0 to 300 mm/Hg
- Measuring accuracy: ± 3 mmHg
- Velcro-closing cuff and inflation bulb, spiral connecting tube
- Wall anchor support with basket bracelet holder
- For professional and non-professional use
- EC medical device

DIR/QM – Rev. 02_02/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.