



OMM2EC

AUTOMATIC SPHYGMOMANOMETER OMRON M2

- Digital Automatic Blood Pressure Monitor
- Compact, fully automatic, user-friendly and quick
- Measures blood pressure and pulse rate
- “IntelliSense” technology for comfortable controlled inflation
- Detects irregular heartbeats
- Memory for 30 sets
- Digital LCD display
- Cuff circumference: 22-32 cm
- Power source: 4 AA batteries 1,5 V (inclusive)
- Contents: monitor, arm cuff, instruction manual, storage case, battery set, blood pressure pass

DIR/QM – Rev. 02_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.