



**TER174**

## **MERCURY-FREE THERMOMETER**

- Clinical mercury-free thermometer for body temperature detection
- The liquid used is a "Gallio/Indio" alloy combination, totally non-toxic and safe for the environment
- Detection range: 35-42° C, precision: 0,10° C
- Dimensions: about 125x12x8 mm
- EC Medical Device
- Measurements must be carried out as part of medical supervision

*DIR/QM - Rev. 02\_04/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.