



FS20C

ADULT FINGER PULSEOXIMETER

- Device for measuring blood oxygen saturation
- Compact and lightweight, equipped with colour LED display
- Detects oxygen saturation (SpO2) and pulse rate
- Battery level indicator
- Provided with 2 AAA batteries
- Dimensions: about 58x33x36mm
- Maximum weight: 60g (with batteries)
- EC medical device
- The measurements must be followed by further specific investigations by competent doctors

DIR/QM – Rev. 02_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.