



HS201C1

HANDHELD SPHYGMOMANOMETER

- Hand-held sphygmomanometer with black cuff
- Large pressure gauge with graduated scale up to 300 mm/Hg
- Nylon cuff with Velcro closing for rapid use
- Inflation bulb with pressure regulation valve
- Measuring accuracy: ± 3 mmHg
- Supplied with transport case
- EC medical device

DIR/QM – Rev. 02_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.