



## HS2000

### SPHYGMOMANOMETER STANDARD BLUE

- Standard aneroid sphygmomanometer
- Pressure gauge with graduated scale up to 300 mm/Hg
- Distinctive colours of the dial: red and black. Red, for immediate recognition of hypo/hypertension states
- Blue nylon cuff with Velcro closing
- Separate inflation bulb with fine pressure regulation valve
- Supplied with transport and protective case
- EC medical device
- Measurements are to be carried out by qualified healthcare professionals

DIR/QM – Rev. 02\_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.