



**HS-201Q1**

## **SPHYGMOMANOMETER PROFESS. I**

- Professional 1 tube sphygmomanometer
- Big manometer with scale up to 300 mm/Hg (70 mm diameter)
- Cuff with Velcro closing
- Deflation valve with screw regulator or rapid deflation button
- Turning manometer that allows the use to both left- and right-handed users
- Practical protection and transport case
- EC medical device
- The device can be used by trained personnel only or doctors

*DIR/QM – Rev. 02\_04/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.