



32565

DUAL HEAD STETHOSCOPE TRAD BLACK

- Traditional line in aluminium, lightweight and latex-free
- Brass headset and anti-shiver chestpiece
- Diaphragm Ø 43,5 mm
- Bell Ø 30 mm
- Tubing: length 62 cm, int. Ø 4 mm, ext Ø 8 mm
- Colour: black
- Class I medical device in accordance with MDR 2017/745

DIR/QM – Rev. 02_02/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.