



**32557**

### **PEDIATRIC STETHOSCOPE TRAD PINK**

- Traditional line in aluminium, lightweight and latex-free. Paediatric
- Brass headset and anti-shiver chestpiece
- Diaphragm Ø 35mm
- Tubing: length 62 cm, int. Ø 4 mm, ext Ø 8 mm
- Colour: pink
- Class I medical device in accordance with MDR 2017/745

*DIR/QM – Rev. 02\_02/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.