



## 32563

## STETHOSCOPE TRAD YELLOW

- Traditional line in aluminium, lightweight and latex-free
- Brass headset and anti-shiver chestpiece
- Diaphragm Ø 43,5 mm
- Tubing: length 62 cm, int. Ø 4 mm, ext Ø 8 mm
- Colour: yellow
- Class I medical device in accordance with MDR 2017/745

DIR/QM - Rev. 02\_02/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.