



**HS-30Q**

## **PROFESSIONAL CHILDREN STETHOSCOPE**

- Professional stethoscope ideal for use in the pediatric sector
- High sensitivity and excellent sound transmission
- The head of the device is made of chromed metal
- Ear tube made of resistant latex-free plastic material, easy to clean and disinfect
- CE medical device in class I according to European Regulation 2017/745 (MDR)
- Imported and distributed by Boscarol Emergency Systems

*DIR/QM – Rev. 03\_18/06/2025*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.