



HS-201X

SPHIGMOMANOMETER PALMAR WHITE

- Palm-sized sphygmomanometer made of white ABS
- Very light medical device that can be used by everyone (about 250 g)
- Large dial with graduated scale up to 300 mm/Hg
- Nylon cuff with Velcro closure for quick application
- Measurement accuracy: ± 3 mmHg on the measurement taken
- Practical protective and transport case
- CE medical device in class Im
- Imported and distributed by Boscarol Emergency Systems

DIR/QM – Rev. 03_18/06/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.