



**3001003**

### **ANEROID SPHYGMO WITH STETHOSCOPE**

- Kit consisting of sphygmomanometer and stethoscope
- Classic sphygmomanometer with nylon cuff, Velcro fastener, latex pear with chrome valve and detachable manometer which can be attached to the cuff
- Chromium-plated air valve with air release control via screw valve
- Scale graduated every 2 mmHg, from 0 to 300 mmHg
- Manual control of air release possible from 260 mmHg to 15 mmHg
- Deflation time max. 10 sec.
- Floating zero (pointer does not have a stopping point but swings freely)
- Supplied with a soft leatherette bag with zip, in a cardboard box
- Colour: black

*DIR/QM – Rev. 01\_12/03/2026*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.