



TR100401

T-TOUCH INFRARED THERMOMETER

- Contact IR technology
- Temperature detected instantly and automatically by placing the thermometer on the forehead.
- Visual fever notification via coloured LED and audible alarm
- Storage of up to 20 results for easy monitoring
- Powered by a non-rechargeable PC 2032 battery (1 piece)
- Automatic shut-off
- Battery life is approximately 1,200 uses
- Comes with protective case
- Dimensions: 2.5 x 4 cm

DIR/QM – Rev. 01_13/03/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.