



**02SF1853EASYICE**

## **NON-WOVEN ICE PACK**

- TNT disposable ice pack for cold therapy with controlled temperature at 0° C for about 30 minutes
- Allows immediate and direct use, preventing swelling and bruising
- Stops small haemorrhages and significantly relieves pain
- Can be stored indefinitely at room temperature
- Pouch made of TNT
- Dimensions: approx. 13.5x17cm
- Pack of 25 pieces

*DIR/QM - Rev. 02\_24/09/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.