



02SF1853EASYICE

NON-WOVEN ICE PACK

- TNT disposable ice pack for cold therapy with controlled temperature at 0° C for about 30 minutes
- Allows immediate and direct use, preventing swelling and bruising
- Stops small haemorrhages and significantly relieves pain
- Can be stored indefinitely at room temperature
- Pouch made of TNT
- Dimensions: approx. 13.5x17cm
- Pack of 25 pieces

DIR/QM - Rev. 02_24/09/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.