



220000066910

## HEIMLICH VALVE KIT

- Non-return valve for thoracic drainage (Heimlich type)
- Valve body in flexible and transparent polyurethane with natural rubber sleeve in "duck's bill" shape acting as a non-return valve
- This device is connected directly onto the patient's pleural or thoracic drain (adult, paediatric), or indirectly, using connector code 800.01 and an extension tube for neonates and infants.
- This product allows thoracic drainage to continue without risk to the patient even when transported, as the non-return valve means that the drain does not have to be clamped.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.