



**250.200**

## **PLEUROTHORAX SET**

- Closed System comprising:
- Secretion Bag 2,0 ltr. with non-return valve, stopcock on the base and three-way stopcock (fixed)
- Transfusion Tube Luer-Lock (unscrewable)
- Puncture Cannula Diameter 2 mm; puncture length 75 mm
- 50-ml-syringe Luer-Lock
- Double and sterile packaging; tear-open package

*DIR/QM – Rev. 02\_25/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.