



39920

POWDER FREE LATEX GLOVES SIZE S 100PCS

- Disposable latex gloves (natural rubber), powder-free with a double internal and external chlorination process (chlorine baths). This process guarantees the absence of powder and reduces considerably protein elements.
- Not sterile, for examination purpose.
- White, versatile ambidextrous with anti-rip wrist protection.
- Micro-textured external surface guarantees high grip and adherence, allowing also maximum tactile sensitivity.
- Suitable for contact with food (oils and fat excluded).
- Use: examination, therapy, diagnostics and laboratory. Aesthetics, cleanliness and hygiene for the community, preparation of drugs and cosmetics, hobby.
- S Size.
- Packaging: 100 pcs.

DIR/QM – Rev. 02_10/10/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.