





42573M

BLACK NITRYL GLOVES SIZE M

- Disposable non-sterile nitrile gloves powder free, for protection, for patients' examination
- Hypoallergenic: latex free, non-allergenic
- Black, ambidextrous with snatch-free wristband and reinforced hems
- The particular formula with which they are made, gives them an exceptional sensibility and softness and at the same time great flexibility. Comparable with latex gloves.
- · Good resistance against chemicals, fuel, oil, grease and detergents
- Micro-coarse surface granting the utmost sensibility
- Food compatible
- Intended use: Best choice for replacing latex gloves in latex-free environment. Ideal for medication, hygiene, laboratory, hobbies and pharmaceutical.
- Size: Medium
- 100 pcs package

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.