



63101

**ALUMINIUM SHEET 60x80 CM**

**ALUMED CLOTH 80x60 cm**

- ALUmed the new sterile medication compress for wounds and burns
- Hypoallergenic, transpiring, antibacterial, absorbent but not wound or burn adherent
- Composed by a layer of aluminium, and with a layer of gauze with a high absorbent power who can provide a considerably reduction of the infection
- The aluminium film, placed directly on the wound or burn, allows the absorption of the sweat or blood without the aluminium sticking on the wound
- It favours epithelization and cicatrisation tissue constitution
- It facilitates medication change thanks to its high absorbing power without adherence
- Dimensions: about 80x60 cm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.