



**CPS258**

## **EMERGENCY EYE WASHING KIT**

- Recommended to all companies and on transport motor vehicles for emergency eye washing and its subsequent medication. Green case with wall support featuring:
- 1 empty eye washing bottle 500 ml.
- 1 eye washing sterile saline solution wash 500 ml
- 2 sterile eye pads
- 2 sterile gauze packages 18x40cm
- 1 pair of lister scissors 14,5cm
- 1 plaster roll 1mx2cm
- 1 towelette package
- 1 hand mirror
- 1 instruction manual
- Dimensions: about 25,5x9x19h cm
- Weight: about 1,2Kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.