



**MAT022**

## **EYE SPLINTER REMOVAL PENCIL**

### **REMOVAL PENCIL**

- Small metallic splinters eye tool remover. One side you have a NEODIMIO magnetic point, which enable you to withdraw all metallic splinters avoiding any eventual dangerous eye approach. On the other side you have an non-abrasive nylon loop which is used for non metallic splinters
- Main body and support made of anodised aluminium and the neoprene seals guaranty a better pencil hold
- Supplied in a practical box
- CE certified
- CND code Q0299 ? RdM code 72774/R

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.