



VAL40950

REIMPLANT LIMBS SET CASE

- Replantation limbs set case, made in shockproof orange ABS, featuring:
- 1 sterile replacement/transport set for finger
- 1 sterile replacement/transport set for hand
- 1 sterile replacement/transport sack for arm
- 1 sterile replacement/transport sack for leg
- 1 haemostatic lace/belt
- 10 sterile gauze packs 10x10cm
- 1 plaster roll 2,5cm
- 1 emergency pad bandage
- 1 Alumed cloth 40x60cm
- 1 anatomic forceps
- 1 surgery forceps
- 1 Pean forceps
- Dimensions: about 42x15x31(h) cm
- Weight: about 7,5 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.