



## **VAL40954**

## **REIMPLANTATION KIT**

- Replantation set in a tearproof impermeable equipped bag with 4 coloured pouches with handles transparent frontal, complete of:
- 1 sterile replacement/transport set for finger
- 1 sterile replacement/transport set for hand
- 1 sterile replacement/transport sack for arm
- 1 sterile replacement/transport sack for leg
- Hemostatic belt/lace
- 1 tourniquet type Martin
- 10 sterile gauze packs 10x10cm
- 1 plaster roll 2,5cm
- 1 emergency pad bandage
- 1 Alumed cloth 40x60cm
- 1 anatomic forceps
- 1 surgery forceps
- 1 Pean forceps
- Dimensions: about 48x16x26(h) cm
- Weight: about 5,5 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.